

Voluson E8

Ultra sound future

Product description

The Voluson® E8 is a premier imaging system delivering a new generation of image quality combined with groundbreaking diagnostic tools and its unique leading volume ultrasound technology.

Highlights

- Matrix array volume technology
- SonoVCAD – Sonography-based Volume Computer Aided Diagnosis
- High resolution transvaginal probe
- High resolution flat panel display
- Electrical height adjustment
- Floating user interface
- TruScan™ architecture
- On-board archive including preview and pre selection



Figure 1. Voluson E8 with Volume Ultrasound technology



General Specifications

Dimensions and Weight

- Height: 1290 mm (50.8 in)
Adjustable: +200 mm (7.9 in)
- Width: 580 mm (22.8 in)
- Depth: 930 mm (36.6 in)
Adjustable: +200 mm (7.9 in)
- Weight (no Peripherals):
120 kg (265 lb)

Electrical Power

- Voltage: 100-130/220-240 Vac
- Frequency: 50/60 Hz
- Power: Max. 1000 VA with on-board Peripherals
- Thermal Output: 3446 BTU/h

Console Design

- 3 Active Probe Ports
(plus 1 non-imaging port)
- Integrated HDD (160 GB)
- Integrated DVD+ R(W) / CD-R(W)
drive
- On-board storage for Peripherals
- Wheels
- Wheel diameter 150 mm
- Integrated locking mechanism that
provides rolling lock
- Integrated cable management
- Front and rear handles

User Interface

Operator Keyboard

- Floating Keyboard:
Rotation: adjustable +/- 40° from
center
Height adjustable + 200mm
- Full-sized, backlit alphanumeric
keyboard
- Ergonomic hard key layout
- Interactive back-lighting
- Integrated recording keys for remote
control of up to 4 Peripherals or
DICOM devices

Touch Screen

- 10.4 in High Resolution color LCD
screen
- Interactive dynamic software menu
- Brightness adjustable

Monitor

- 15" LCD monitor
- High brightness with 350 cd/m²
- Tilt/Rotate Adjustable Monitor
Tilt Angle: + 10°/- 90°
Rotate Angle: 360°
- Digital brightness & contrast
adjustment

System Overview

Applications

- Abdominal
- Obstetrical
- Gynecological
- Small parts
- Vascular / Peripheral
- Pediatric and Neonatal
- Urological
- Cardiology
- Neurology
- Orthopedic

Operating Modes

- B-Mode (2D)
- M-Mode (M)
- M-Color-Mode(MC)
- Color Flow Mode (C)
- Power Doppler Imaging (PD)
- Tissue Doppler Imaging (TD)
- HD-Flow Imaging (HD-Flow)
- PW Doppler with high PRF (PW)
- B-Flow (BF)
- Extended View (XTD View)
- Coded Contrast Imaging
(Contrast Media)
- Volume Mode (3D/4D):
 - 3D Static
 - 4D Real Time
 - VCI-A, VCI-C
 - STIC /Color, Angio, HD-Flow,
Contrast & B-Flow
 - 4D Biopsy

Scanning Methods

- Electronic Sector
- Electronic Convex
- Electronic Linear
- Mechanic Volume Sweep

Transducer Types

- Sector Phased Array
- Convex Array
- Microconvex Array
- Linear Array
- Volume probes `4D`:
 - Convex Array
 - Microconvex Array
 - Linear Array
 - 1.5D Convex array
 - 1.5D Linear array
 - 1.25D Convex Array
 - 1.25D Linear Array

System Standard Features

- State-of-the-art user interface with
high resolution 10.4 inch LCD touch
panel
- Automatic Tissue Optimization
- Tissue Doppler
- Coded Harmonic Imaging
- Coded Excitation (CE)
- HD-Flow
- XTD
- SRI III (Speckle reduction imaging)
- CrossXBeamCRI (Compound
Resolution Imaging)
- Static 3D Mode:
 - B Mode only

- B + Power Doppler Mode
- B + CFM Doppler Mode
- B + HD-Flow Mode
- B + CRI
- B + CRI + CFM
- B + CRI + PD
- B + CRI + HD-Flow
- B + Contrast
- B + B-Flow
- Focus&Frequency Composite (FFC)
- High Resolution Zoom
- Pan Zoom
- Steering
- Virtual Convex
- Beta-View
- Patient information database
- Image Archive on hard drive
- 3D/4D data compression
(lossy/lossless)
- Inversion
- Real-time automatic Doppler calcs
- Measurement & Calculations
including Worksheets/Report for :
 - OB
 - GYN
 - Vascular
 - Cardio
 - Abdominal
 - Small-Parts
 - Urology
 - Pediatrics
 - Ortho
 - Neurology
- Multigestational Calculations

System Options

- 4D Real Time
- VOCAL II
- VCI (Volume Contrast Imaging)
- SonoVCAD
- DICOM
- 4D – STIC:
 - STIC
 - STIC + Power Doppler Mode
 - STIC + CFM Doppler Mode
 - STIC + HD-Flow Mode
 - STIC + CRI
 - STIC + CRI + CFM
 - STIC + CRI + PD
 - STIC + CRI + HD-Flow
 - STIC + Contrast
 - STIC + B-Flow
- B-Flow
- T.U.I – Tomographic Ultrasound
Imaging
- Coded Contrast Imaging
- Foot Switch, with programmable
functionality

Peripheral Options

- Integrated printers:
 - B&W thermal printer
 - Color thermal printer
- DVD Recorder
- ECG Digital Module
- External Color PC desktop printer &
connection kits

Display Modes

- Simultaneous Capability
 - B+PW
 - B+CFM, B+PD, B+TD, B+HD-Flow
 - B+M
 - B+3D, B+4D
 - B+CRI
 - B+SRI
 - B+CRI+SRI
 - BF+SRI
 - Contrast+SRI
 - B+CRI/3D+CRI
 - B+SRI/3D+SRI
 - B+CRI+SRI/3D+CRI+SRI
 - B+CRI/4D+CRI
 - B+SRI/4D+SRI
 - B+CRI+SRI/4D+CRI+SRI
 - B+CRI/STIC+CRI
 - B+SRI/STIC+SRI
 - B+CRI+SRI/STIC+CRI+SRI
 - B/B+CRI
 - B/B+SRI
 - B/B+SRI+CRI
 - B/CFM+CRI
 - B/CFM+SRI
 - B/CFM+CRI+SRI
 - B/PD+CRI
 - B/PD+SRI
 - B/PD+CRI+SRI
 - B/HD-Flow + CRI
 - B/HD-Flow +SRI
 - B/HD-Flow +CRI+SRI
- Real-time Triplex Mode
 - B/CFM/PW
 - B/PD/PW
 - B/HD-Flow/PW
- Selectable alternating Modes
 - B+PW
 - B/CFM+PW
 - B/PD+PW
 - B/HD-Flow+PW
 - B+CFM or PD or HD-Flow
- Multi-image (split, quad)
 - Live and/or frozen
 - split: B+B, B/CFM + B/CFM, or B/PD or B/TD or B/HD-Flow or BF+BF, Contrast + Contrast
 - split: B+B/CFM or PD or HD-Flow
 - split: B+PW or M
 - split: Frame Review / XTD-View
 - quad: B+B+B+B or BF or Contrast, B/CFM+B/CFM+B/CFM +B/CFM or B/PD or B/TD or B/HD-Flow
 - Independent Cine playback
 - Quad: A+B+C+3D or 4D
 - 3x3: T.U.I Overview + 8 parallel slices
 - Quad: T.U.I Overview + 3 parallel slices
 - Split: T.U.I Overview + 1 slice
- Zoom Read/Write (with or without overview image)
- Colorized Image
 - Colorized B
 - Colorized M
 - Colorized PW
 - Colorized 3D
- Time line display
 - Independent Dual B/PW Display

- Display Formats
 - Top/ Bottom selectable format (Size: 1/2:1/2; 1/3:2/3; 2/3:1/3)
 - Left/Right selectable format (Size: 1/2:1/2; 1/3:2/3; 2/3:1/3)

Display Annotation

- Patient Name:
 - Last: max 32 characters
 - First: max 15 characters
 - Middle: max 15 characters
- ID: max 32 characters
- Accession #: max 16 characters
- Hospital Name: max 30 Characters
- Sonographer (up to 5 characters are displayed depending on font size)
- Gestational age (OB) or LMP (Gyn)
- Birth date (selectable)
- Date: 3 Types selectable
 - MM/DD/YYYY
 - DD/MM/YYYY
 - YYYY/MM/DD
- Time: 2 types selectable
 - 24 hours
 - 12 hours
- Probe Name
- Application Name
- Gray Scale bar
- Depth Scale
- Focal Zone Marker
- Frame Rate
- Zoom Start/Depth
- B-Mode
 - User program
 - Receiver Frequency
 - Acoustic Power
 - Gain
 - Dynamic Contrast
 - Gray Map
 - Edge Enhance
 - Persistence
 - SRI, CRI
 - Focal Zone Markers
 - Depth Scale Marker
 - Probe Orientation
- M-Mode
 - Gain
 - Dynamic Contrast
 - Edge Enhance
 - Reject
 - M-Cursor
 - Time Scale
- Doppler Mode
 - Acoustic Power
 - Gain
 - Angle
 - Sample Volume Depth and Width
 - Wall Filter
 - Velocity or Frequency Scale
 - Spectrum Inversion
 - Time Scale
 - PRF
 - HPRF
 - Doppler Frequency
- Color Flow Imaging Modes (CFM, PD, TD, HD-Flow)
 - Acoustic Power
 - Color Gain
 - Color Balance

- Color Balance Marker
- Quality
- Wall Motion Filter
- PRF
- Color Map
- Color Scale: 2 types
 - Power and Symmetrical Velocity Imaging
- Color Velocity Range
- Spectrum Inversion
- 3D/4D Mode
 - 3D/4D Sub Program
 - Threshold
 - Quality
 - Volume Box Angle
 - Mix
 - Acquisition Mode
 - Compression
 - Orientation Markers
 - T.U.I.: slice distance (0.5-10mm)
 - T.U.I.: slice position in overview image
 - SonoVCAD
- TGC Curve
- Cine Frame Number
- Recorder Status
- Body Pattern: 111 types organized in 10 anatomical groups
- Measurement Results
- Displayed Acoustic Output
 - TIS: Thermal Index Soft Tissue
 - TIC: Thermal Index Cranial (Bone)
 - TIB: Thermal Index Bone
 - MI: Mechanical Index
- Power output (%)
- Biopsy Guide Line
- ECG Line
- Trackball function (Trackball and Trackball buttons)
- GE Logo
- Zoom overview image (zoom box position)

System Parameters

System Setup

- Pre-programmable Categories date format.
- User Programmable Preset Capability, User program etc.
- Languages: English, French, German, Spanish, Italian + additional languages loadable
- EUM Languages: English, German, Spanish, Portuguese, Italian, French, Mandarin Chinese, Japanese
- Up to 400 Programmable Annotations organized in 10 anatomical groups

Measure Setup

- M&A Setup including Add, Delete, Edit and Reorder of measure items
- Application Setup including several parameters of Measurement, Doppler Trace and Calculation presets

- Global Setup including several parameters of Measurement, Cursor and Result window presets

Pre-Processing

- Write Zoom up to 8x
- B/M-Mode
 - Gain
 - TGC
 - Dynamic Range
 - Acoustic Output
 - Transmission Focus Position
 - Transmission Focus Number
 - Transmission Frequency
 - Edge Enhancement
 - Persistence Control
 - Line Density Control
 - Reject
 - Sweep Speed
 - M-Cursor position
- PW-Mode
 - Gain
 - Dynamic Range
 - Acoustic Output
 - Transmission Frequency
 - PRF
 - Wall Filter
 - Sample Volume Gate Length, Depth, Pos
 - Velocity Scale
 - Sweep Speed
- Color Flow Imaging Modes (CFM, PD, TD, HD-Flow)
 - Gain
 - Acoustic Output
 - PRF
 - Wall Motion Filter
 - Line density
 - Ensemble
 - Dynamic
 - Smooth (Rise and Fall)
 - Frequency
 - Balance
 - Line Filter
 - Quality
 - Artifact Suppression

Post-Processing

- Read Zoom: 0.8x - 3.4x Zoom (with HD-Zoom functionality up to 22x Zoom)
- B/M-Mode
 - Gray Map
 - Colorized B and M
 - Speckle Reduction Imaging (SRI III)
- PW Mode
 - Gray Map
 - Baseline Shift
 - Angle Correction
 - Colorized D
 - Scale (KHz, m/s, cm/s)
 - Trace
 - Invert
- Color Flow Imaging Modes (CFM, PD, TD, HD-Flow)
 - Color Map
 - Display Threshold
 - Display Mode (V, V-T,T,P,P-T) (CFM only)

- Scale (CFM and HD-Flow)
- Baseline
- B-Flow
- Gray Map

Image Processing and Presentation

- Digital Beamformer
- 16896 system processing channel technology
- Displayed Imaging Depth: 0 – 30 cm
- Minimum Depth of Field: 0 – 1 cm (Zoom, probe dependent)
- Maximum Depth of Field: 0 – 30 cm (probe dependent)
- Transmission Focus
 - 1- 5 Focus Points selectable (probe and application dependent)
 - Focal Zone position, up to 7 steps
- Continuous Dynamic Receive Focus / Continuous Dynamic Receive Aperture
- 256 shades of gray
- 16,8 Mio Colors 24 bit
- Up to 180 dB Dynamic. Range adjustable by selecting 12 Dynamic Contrast Curves
- Image Reverse: Right/ Left
- Rotation: 0°, 180°

CINE Memory/Image Memory

- CINE Memory: up to 256 MB (up to 3000 2D images)
- Dual Image CINE Display
- Quad Image CINE Display
- CINE image number display
- CINE Review Loop
- CINE Review Speed: 4 speeds: 25/50/100/200%
- Length of CINE Sequence Review selectable (start/end image)
- Measurements/ Calculations & Annotations on CINE Playback

Image/Volume Storage

- On-board data storage software:
 - Image data stored as:
 - Raw Data file (proprietary format) or
 - DICOM file (Single- or Multiframe)
 - Volume file:
 - Format: proprietary
 - Size: typically: 0.8 - 5MB (depending on probe and adjusted volume size)
 - Lossy and lossless compression available.
 - Typical compression rates are 50% with lossless compression, 15% with lossy compression but maximum quality and 5% with lossy compression and reduced quality (approximate values).
 - Cine Review

- Single Volume (raw data, conversion to Cartesian format)
- Volume Cine (raw data)
- 2D Raw Data (incl. Doppler and Motion Data)
- 3D Movie
- Measure Reports
- Information from past exams
- Export functions:
 - Format:
 - BMP, TIFF or JPEG
 - Sequence of BMPs
 - AVI or Quicktime
 - Voluson Format (proprietary)
 - DICOM Format
 - Volumes/Raw Data
 - Export to: DVD+ R(W) / CD-R(W), Network, USB devices
- Backup function to DVD+ R(W) / CD-R(W), Network, USB devices
- AVI-Files and Quicktime-Files: conversion and export to: DVD+ R(W) / CD-R(W), Network, USB devices
- Compression available
- Hard Drive Data Storage: 100 GB
- Compare old images with current exam
- Reload of archived data sets

Connectivity

- Ethernet network connection
- USB for USB devices
- DICOM support (option)
 - Verify
 - Print
 - Store
 - Modality Worklist
 - Structured Reporting
 - Storage Commitment
 - MPPS (Modality performed procedure step)
 - Media Exchange
 - Off network / mobile storage queue

Scanning Parameters

B-Mode

- B Acoustic Power: 1-100%
- B Gain: +/-15dB range, 1dB steps
- Slide pots: +/- 15dB
- Dynamic Range: max.180 dB, 12 dynamic Contrast curves
- Persistence: 8 steps
- B Gray Scale Map: 12 maps
- B Edge Enhancement: 5 steps
- Line Filter: 3 steps
- Reject: range 0-255, step size 5,
- Frequency Selection: 3 steps (multi-frequency, wideband probes)
- Quality (Line Density): 3 steps
- Scanning Size (FOV or Angle Depending on probe)
- B Colorization: 6 chroma maps
- BetaView (Volume probes only)
- 3 Display Formats (60:40, 50:50, 40:60)

M-Mode

- M Acoustic Power: 1-100%
- M Gain: +/-15dB range, 1dB steps
- Slide pots: +/- 15dB
- Dynamic Range: max.180 dB, 12 Dynamic Contrast curves
- M Gray Scale Map: 12 maps
- M Edge Enhancement: 5 steps
- M Sweep Speed: 4 types
- M Colorization: 6 chroma maps
- M Reject: range 0-255, step size 5,

M-Color Flow Mode

- Acoustic MCFM Power: 1-100%
- Frequency range: 1-15Mhz (Depending on the probe, 3 steps high, mid, low)
- MCFM Color Maps: 8 maps
- CFM Gain: +/-15dB range, 1 dB steps
- CFM Velocity Scale Range: PRF: 100Hz to 13kHz
- Wall Filter: 8 – 3000 Hz
- Ensemble (color shots per line) 8-16, step size 1
- Gentle color filter
- Smooth filter:
 - Rise: 12 steps
 - Fall: 12 steps
- CFM Spectrum Inversion
- CFM Baseline Shift: 17 steps
- Pre-settable and independently adjustable B-, M and MCFM Gain
- CFM Threshold: 1 – 255 steps
- Balance: 25 – 225, step size 5
- Artifact suppression: on/off
- Color Display Mode:
 - V (Velocity)
 - V-T (Velocity + Turbulence)
 - V-P (Velocity + Power)
 - T (Turbulence)
 - P-T (Power + Turbulence)
- Real-time Triplex Mode: B + M +MCFM in any depth

Spectral Doppler Mode (PW)

- Acoustic Power: 1-100%
- Transmit Frequency Range:
 - PW: 1 – 15Mhz
- Gain: +15/-25dB range, 1dB steps
- Displayed Dynamic Range: 10 – 40 dB, 2 dB steps
- Gray Scale Map: 12 maps
- PW Wall Filter: 70 – 500Hz, 7 steps, PRF dependent
- Colorization: 6 chroma maps
- PW PRF: 1.3 – 22.0 kHz
- PW: Velocity Scale Range (Depending on the probe Frequency)
 - 2MHz, 0°, max. zero shift range: 1cm/s - 8m/s
 - 2MHz, 60°, max. zero shift range: 1cm/s - 16m/s
- PW Sweep Speed:
 - Time Resolution: Simplex 2.2, 3.3, 4.4, 6.6,10 msec Duplex/Triplex 4.4, 6.6,10 msec.
- Sample Volume Length: (.7mm, 1-10 (steps 1mm), 15mm)

- Spectrum Analyzer (FFT):
 - max: 256 channels
 - 255 amplitude levels
- Angle Correction: ± 0-85°, 1° step Available before Freeze and after Freeze
- Steered Linear: 0° - 25° (Depending on probe)
- Spectrum Inversion
- Baseline Shift: +/- 8 steps from center
- Doppler Auto Trace

Color Flow Mode

- Acoustic Power: 1-100%
- Frequency range: 1-15Mhz (Depending on the probe, 3 steps high, mid, low)
- Color Map: 8 maps
- CFM Gain: +/-15dB range, 0.2 dB steps
- CFM Velocity Scale Range: PRF: 100Hz to 11kHz
 - less than +/- 0.3 cm/s
 - max: +/- 5.5 m/s
- Wall Filter: 8 – 3000 Hz, 7 steps
- Ensemble (color shots per line) 7-31, step size 1
- Line Density: 10 steps
- Gentle color filter
- Line Filter: 8 steps
- Smooth filter:
 - Rise: 12 steps
 - Fall: 12 steps
- CFM Window Size:
 - max: same as B-image size
 - +/- 25 ° (probe dependent)
- CFM Spectrum Inversion
- CFM Baseline Shift: 17 steps
- Pre-settable and independently adjustable B-Mode Gain in B/CFM-Mode
- CFM Threshold: 1 – 255 steps
- Balance: 25 – 225, step size 5
- Artifact suppression: on/off
- Color Display Mode:
 - V (Velocity)
 - V-T (Velocity + Turbulence)
 - V-P (Velocity + Power)
 - T (Turbulence)
 - P-T (Power + Turbulence)
- Real-time Triplex Mode: B + CFM/PW in any depth
- 3 Scales (kHz, cm/s, m/s)

Power Doppler Imaging (PD)

- Acoustic Power: 1-100%
- Frequency range: 1-15Mhz (Depending on the probe, 3 steps high, mid, low)
- PD Color Maps: 8 maps
- Gain: +/-15dB range, 0.2 dB steps
- Velocity Scale Range:
 - PRF: 100Hz to 11kHz
- Wall Filter: 8 – 3000 Hz, 7 steps
- Ensemble (color shots per line) 7-31, step size 1
- Line Density: 10 steps

- Gentle color filter
- Line Filter: 8 steps
- Smooth filter:
- Rise: 12 steps
- Fall: 12 steps
- PD Window size:
 - max.: same as B-image size
- Maximum Steer-able Angle
 - +/- 25 ° (probe dependent)
- Pre-settable and independently adjustable B-Mode Gain in B/PD-Mode
- PD Threshold: 0 – 255 steps
- Artifact suppression: on/off
- Balance: 25 – 225, step size 5
- Real-time Triplex Mode: B + PD/PW in any depth

HD-Flow

- Acoustic Power: 1-100%
- Frequency range: 1-15Mhz (Depending on the probe, 3 steps high, mid, low)
- HD Color Map: 8 maps
- HD Gain: +/-15dB range, 0.2 dB steps
- HD Velocity Scale Range: PRF: 100Hz to 11kHz
 - less than +/- 0.3 cm/s
 - max: +/- 5.5 m/s
- Wall Filter: 8 – 3000 Hz, 7 steps
- Ensemble (color shots per line) 7-31, step size 1
- Line Density: 10 steps
- Line Filter: 8 steps
- Smooth filter:
 - Rise: 12 steps
 - Fall: 12 steps
- HD Window Size:
 - max: same as B-image size
 - +/- 25 ° (probe dependent)
- HD Spectrum Inversion
- HD Baseline Shift: 17 steps
- Pre-settable and independently adjustable B-Mode Gain in B/HD-Mode
- HD Threshold: 1 – 255 steps
- Balance: 25 – 225, step size 5
- Artifact suppression: on/off
- Real-time Triplex Mode: B + HD-Flow/PW in any depth

Tissue Doppler Imaging (TD)

- Acoustic Power: 1-100%
- Frequency range: 1-15Mhz (Depending on probe, 3 steps high, mid, low)
- TD Color Maps: 4 maps
- Gain: +/-15dB range, 1 dB steps
- Velocity Scale Range: PRF: 100Hz to 11kHz
- Wall Filter: 8 – 3000 Hz, 7 steps
- Ensemble (color shots per line) 7-31, step size 1
- Line Density: 10 steps

- Gentle color filter
- Line Filter: 8 steps
- Smooth filter:
- Rise: 12 steps
- Fall: 12 steps
- TD Window size:
- Max: same as B-image size
- CFM Spectrum Inversion
- Pre-settable and independently adjustable B-Mode Gain in B/TD-Mode
- TD Threshold: 1 – 255 steps
- Balance: 25 – 225, step size 5
- Real-time Triplex Mode: B + TD/PW in any depth

Auto Optimization

- Available in:
 - B-Mode
 - PW Doppler

Coded Excitation (CE)

- Available on the following probes:
 - M6C-D
 - RAB4-8P
 - RIC6-12-D
 - RIC5-9-D

Coded Harmonic Imaging

- Harmonic Imaging
- Available on all probes

Compound Resolution Imaging (CRI)

- CRI
- 1-8 steps selectable
- Available on all probes

Focus Frequency Composite (FFC)

- Available on the following probes:
 - M6C-D
 - 4C-D
 - RAB2-5-D
 - RAB4-8-D
 - RIC6-12-D
 - RIC5-9-D

Speckle Reduction Imaging (SRI III)

- 1-12 steps selectable
- Available on all probes

Volume Mode (3D/4D)

- Acquisition Modes:
 - 3D Static: B-Mode (incl. CRI)
 - 3D B-Flow (optional)
 - 3D Angio: B/Power Doppler (incl. CRI)
 - 3D CFM: B/Color Doppler (incl. CRI)
 - 3D HD-Flow: B/HD-Flow (incl. CRI)
 - 3D Contrast: B/Contrast
 - 4D Real Time: (optional)
 - 4D Biopsy (optional)
 - VCI-A, VCI-C (optional)
 - STIC: Fetal Cardio (optional) (incl. CRI)
 - STIC Angio: B/Power Doppler

- (optional) (incl. CRI)
- STIC CFM: B/Color Doppler (optional) (incl. CRI)
- STIC HD-Flow: B/HD-Flow (optional) (incl. CRI)
- STIC B-Flow (optional)
- Visualization Modes:
 - 3D Rendering (diverse surface and intensity projection modes)
 - Sectional Planes (3 Section planes perpendicular to each other)
 - Niche: 3DStatic only
 - VOCAL II (optional): semi-auto/manual segmentation tool (segmentation using touch screen), (3D Static only) + Threshold Volume: measure volume below and above a threshold
 - T.U.I.(optional): Tomographic Ultrasound Imaging (overview image + parallel slices)
 - VCI Static (optional): 3D Static only (Sectional planes with VCI)
 - SonoVCAD Heart (optional): Sonography-based Volume Computer Aided Diagnosis
- Render Mode:
 - Surface texture, Surface Smooth, max-, min- and X-ray (average intensity projection), Gradient, Inversion, Glass Body, Mix Mode of two render Modes
- 3D Movie
- Curved render start
- SRI post-processing for A, B and C-Plane and rendered image
- MagiCut: 3D/4D Cut tool
- Display Format:
 - Quad: A-/B-/C-Plane/3D
 - Dual: A-Plane/3D
 - Dual: A-Plane/VCI
 - Single: 3D or A- or B- or C-Plane
 - TUI 3x3: Overview image + 8 slices
 - TUI Quad: Overview image + 3 slices
 - TUI Dual: Overview image + 1 slice
 - TUI Single: slice
- 4D Volume Frames/sec max: 40
- 4D Volume Cine: 128 Volumes

Virtual Convex

- Provides a convex field of view for all linear transducers
 - SP10-16-D
 - RSP6-16-D

Measurements / Calculations

Generic B-Mode and 3D

- Distance
 - Distance (Point to Point)
 - Distance (Line to Line)

- 2D Trace (Trace Length)
- Stenosis (% Dist)
- 2D Trace (Point Length)
- Area/Circumference
 - Ellipse
 - Trace (Line)
 - Trace (Point)
 - Stenosis (% Area)
 - Area 2 Dist.
- Volume: following Methods:
 - 1 Distance
 - 1 Ellipse
 - 1 Dist. + Ellipse
 - 3 Distance
 - Planimetric Volume (3D only)
- Angle:
 - Angle (3 Point)
 - Angle (2 Line)

Generic M-Mode

- Distance (Point to Point)
- Time
- Slope
- HR
- Stenosis (% Dist)

Generic Doppler

Measurements / Calculations

- Auto & Manual Trace:
 - PS (Peak Systole)
 - ED (End Diastole)
 - MD (Min. Diastole)
 - PS/ED (Ratio)
 - PI (Pulsatility Index)
 - RI (Resistance Index)
 - TAmx (Time avg. max. Velocity)
 - TAmn (Time avg. mean Velocity)
 - VTI (Velocity Time Integral)
 - Heart Rate
- Single Measurements:
 - Velocity
 - Acceleration,
 - RI,
 - PI,
 - PS/ED,
 - Time,
 - HR
- PG (Pressure Gradient)
 - PG mean
 - PG max

Real-time Doppler Auto

Measurements / Calculations

- PS (Peak Systole)
- ED (End Diastole)
- PD (Peak Diastole)
- MnV (Mean Velocity)
- VTI (Velocity Time Integral)
- RI (Resistance Index)
- PI (Pulsatility Index)
- S/D (Ratio)
- HR (Heart Rate)

Obstetrics

Measurements / Calculations

- Gestational Age by:
 - AC (Abdominal Circumference)

- APTD (Anterior Posterior Thoracic Diameter)
 - APTDxTTD
 - BOD (Binocular Distance)
 - BPD (Biparietal Diameter)
 - CEREB (Cerebellum)
 - CLAV (Clavicula)
 - CRL (Crown Rump Length)
 - EFW (Estimated Fetal Weight)
 - FIB (Fibula)
 - FL (Femur Length)
 - FTA (Fetal Trunk Area)
 - GS (Gestational Sac)
 - HC (Head Circumference)
 - HL (Humerus Length)
 - LV (Length of Vertebra)
 - MAD (Middle Abdomen Diameter)
 - OFD (Occipital Frontal Diameter)
 - RAD (Radius)
 - TAD (Transverse Abdominal Diameter)
 - TCD (Transverse Cerebellar Diameter)
 - TIB (Tibia Length)
 - TTD (Transverse Thoracic Diameter)
 - ULNA (Ulna Length)
 - UmbArt RI, PI (Umbilical Artery)
 - UtArt RI, PI (Uterine Artery)
 - Use measurement results from other systems for fetal trending (past exam)
 - Estimated Fetal Weight (EFW) by:
 - AC
 - AC, BPD
 - AC, BPD, FL
 - AC, BPD, FL, HC
 - AC, FL
 - AC, FL, HC
 - BPD, FTA, FL
 - BPD, MAD, FL
 - BPD, TTD
 - BPD, APTD, TTD, FL
 - BPD, APTD, TTD, LV
 - Calculations and Ratios
 - FL/BPD
 - FL/AC
 - FL/HC
 - HC/AC
 - CC/TC
 - CI (BPD/OFD) (Cephalic Index)
 - Va/Hem, Vp/Hem
 - AFI (Amniotic Fluid Index)
 - Tables / Calculations by:
 - ASUM, Bahlman, Baschat, Brenner, Bunduki, Campbell, Crequat, Chitty, Daya, Doubilet, Eik-Nes, Goldstein, Hadlock, Hansmann, Hellman, Hill, Hobbins, Hohler, Holländer, Jeanty, Johnsen, JSUM, Kurmanavicius, Kurtz, Lessoway, Mari, Marsal, Merz, Nelson, Nicolaides, O'Brien, Okai, Osaka, Persson, Persutte, Rempen, Robinson, Sabbagha, Schild, Shephard, Shinozuka, Sonek, Tokyo University, Warda, Williams, Yarkoni
 - Programmable OB Tables
 - Programmable OB Formuals
 - OB Report including:
 - Measure results (Calc)
 - Measure results (Generic)
 - Fetal Qualitative Description (Anatomical survey)
 - Fetal Environmental Description (Biophysical profile)
 - Fetal Graphical Trending
 - Fetal Compare
- ## Gynecology
- ### Measurements / Calculations
- Right/Left Ovary Length, Width, Height, Volume
 - Right/Left Kidney Length, Width, Height, Volume
 - Uterus Length, Width, Height, Volume
 - ENDO (Endometrial thickness)
 - Cervix Length
 - Follicular measurements (12)
 - Ovarian Artery
 - Uterine Artery
 - Vessel
 - Summary Reports
 - Possibility to add OB measurements to Gyn application
- ## Vascular
- ### Measurements/Calculations
- CCA (Common Carotid Artery)
 - ICA (Internal Carotid Artery)
 - ECA (External Carotid Artery)
 - Vertebral Artery
 - Subclavia
 - Bulb
 - Vessels
 - Summary Reports
- ## Neurology
- ### Measurements/Calculations
- ACA (Anterior Cerebral Artery)
 - MCA (Middle Cerebral Artery)
 - PCA (Posterior Cerebral Artery)
 - Basilar Artery
 - A-Com. A (Anterior Common Artery)
 - P-Com. A (Posterior Common Artery)
 - CCA (Common Carotid Artery)
 - ICA (Internal Carotid Artery)
 - Vertebral Artery
 - Vessels
 - Summary Reports
- ## Cardiology
- ### Measurements / Calculations
- 2D Mode:
 - Simpson (Single & Bi-Plane)
 - Volume (Area Length)
 - LV-Mass (Epi & Endo Area, LV Length)
 - LV (RVD, IVS, LVD, LVPW)
 - LVOT Diameter
 - RVOT Diameter
 - MV (Dist A, B, Area, PISA)
 - TV (Diameter)
 - AV/LA (Ao & LA Diam.)
 - PV (Diameter)
 - M Mode:
 - LV (IVS, LVD, LVPW, RVD)
 - AV/LA (Ao Diam, LA Diam, AV Sep., AoRoot Ampl.)
 - MV (D-E, E-F Slope, A-C Intervall, E-EPSS, E-S Dist.)
 - HR (Heart Rate)
 - Spectral Doppler Mode:
 - MV (Mitral Valve)
 - AV (Aortic Valve)
 - TV (Tricuspid Valve)

- PV (Pulmonary Valve)
- LVOT & RVOT-Doppler (Left & Right Ventricle Outflow Tract)
- Pulmonic Veins
- PAP (Pulmonary Artery Pressure measurement)
- HR (Heart Rate)
- R-R Interval
- CFM Mode:
 - PISA-Radius
 - PISA-Alias Velocity

- Additional Calculations
 - Diast. Vol.(Bi)
 - Syst.Vol.(Bi)
 - Stroke Volume
 - Cardiac Output
 - Eject. Fraction
 - Fract. Shortening FS
 - Myocardial Thickness
 - LA/Ao, Ratio, E/A
 - Peak Gradient, Peak Gradient Acceleration
 - Mean Gradient, Mean Gradient Acceleration
 - VTI, TVA, PG, PHT, MVA, AVA, ERO, etc

Abdominal

Measurements/Calculations

- 2D Mode:
 - Liver
 - Gallbladder
 - Pancreas
 - Spleen
 - Left/Right Kidney
 - Renal Artery
 - Aorta
 - Portal Vein
 - Vessel
- M Mode:
 - Renal Artery
 - Aorta
 - Vessel
- Spectral Doppler Mode:
 - Renal Artery
 - Aorta
 - Vessel
 - Portal Vein
- Summary Reports

Small Parts

Measurements/Calculations

- Thyroid
- Testicle
- Vessel
- Summary Reports

Urology

Measurements/Calculations

- Bladder
- Prostate
- Testicle
- Left/Right Kidney
- Renal Artery
- Vessel

- Summary Reports including PSAD, PPSA(1), PPSA(2) calculation

Pediatric

Measurements/Calculations

- Hip Joint

Probes

- 4C-D Wide Band Convex Probe
 - Applications: Abdomen, OB, GYN
 - Maximum Band Width (-20dB): 1.5 - 4.6 MHz
 - Number of Elements: 128
 - Convex Radius: 60 mm
 - FOV: 58°
 - Foot Print: 60.8 x 13 mm
 - Doppler Transmission
 - Low Freq.: 2,00 MHz
 - Mid Freq.: 2,92 MHz
 - High Freq.: 3,33 MHz
 - Harmonic Transm. Frequency:
 - Resolution: 2.32 MHz
 - Normal: 2.22 MHz
 - Penetration: 2.00 MHz
 - Biopsy Guide available: 4C, Multi-Angle, disposable with reusable bracket
- M6C Wide Band Convex Probe (1.25D Array)
 - Applications: Abdomen, OB, GYN, Pediatrics
 - Maximum Band Width (-20dB): 2.14-6.10 MHz
 - Number of Elements per row: 192
 - Number of rows: 5
 - Convex Radius: 50 mm
 - FOV: 60°
 - Foot Print: 55x18mm
 - Doppler Transmission
 - Low Freq.: 2,78 MHz
 - Mid Freq.: 3,45 MHz
 - High Freq.: 4,01 MHz
 - Harmonic Transm. Frequency:
 - Resolution: 3,04 MHz
 - Normal: 2,78 MHz
 - Penetration: 2,44 MHz
 - Biopsy Guide available: M7C, Multi-Angle, disposable with reusable bracket
- SP10-16-D Wide Band Linear Probe
 - Applications: Small Parts, Peripherals. Vascular, Pediatrics, Ortho
 - Maximum Band Width (-20dB): 4.5 - 16.5 MHz
 - Number of Elements: 192
 - FOV: 33.7 mm
 - Foot Print: 38 x 5.0 mm
 - Doppler Transmission
 - Low Freq.: 7,70 MHz
 - Mid Freq.: 8,34 MHz
 - High Freq.: 9,10 MHz
 - Harmonic Transm. Frequency:
 - Resolution: 6,45 MHz
 - Normal: 5,70 MHz
- RAB2-5-D Wide Band Convex Volume Probe
 - Applications: Abdomen, OB, GYN
 - Maximum Band Width (-20dB): 2 - 5 MHz
 - Number of Elements: 192
 - Convex Radius: 40,5 mm
 - Volume Sweep Radius: 20,15mm
 - FOV: 80° (B), 85° x 80° (Volume scan)
 - Foot Print: 53,2 x 40,6
 - Doppler Transmission
 - Low Freq.: 2,38 MHz
 - Mid Freq.: 2,95 MHz
 - High Freq.: 3,33 MHz
 - Harmonic Transm. Frequency:
 - Resolution: 2,13 MHz
 - Normal: 2,00 MHz
 - Penetration: 2,00 MHz
 - Biopsy Guide Available: PEC74, Single-Angle, Reusable
- RAB4-8-D Wide Band Convex Volume Probe
 - Applications: Abdomen, OB, GYN, Pediatric, Urology
 - Maximum Band Width (-20dB): 4 - 8,5 MHz
 - Number of Elements: 192
 - Convex Radius: 41.6 mm
 - Volume Sweep Radius: 19,95 mm
 - FOV: 70° (B), 85° x 70° (Volume scan)
 - Foot Print: 53,2 x 40,6 mm
 - Doppler Transmission
 - Low Freq.: 3,04 MHz
 - Mid Freq.: 3,45 MHz
 - High Freq.: 4,01 MHz
 - Harmonic Transm. Frequency:
 - Resolution: 3,23 MHz
 - Normal: 2,78 MHz
 - Penetration: 2,57 MHz
 - Biopsy Guide Available: PEC74, Single-Angle, Reusable
- RIC5-9-D Wide Band Convex Volume Probe
 - Applications: OB, GYN, Urology
 - Band Width (-20dB): 3.7-9.3MHz
 - Number of Elements: 192
 - Convex Radius: 11.6 mm
 - Volume Sweep Radius: 11.6 mm
 - FOV: 146° (B), 146°*90° (Volume scan)
 - Foot Print: 32 x 27 mm
 - Doppler Transmission
 - Low Freq.: 5,26 MHz
 - Mid Freq.: 5,88 MHz
 - High Freq.: 6,66 MHz
 - Harmonic Transm. Frequency:
 - Resolution: 4,30 MHz
 - Normal: 4,30 MHz
 - Penetration: 3,84 MHz
 - Biopsy Guide Available: PEC63, Single-Angle, Reusable

- Penetration: 5,27 MHz
- Steered Angle: Max. 25°
- Biopsy Guide available: PEC64

- RIC6-12-D Wide Band Convex Volume Probe
 - Applications: OB, GYN, Urology
 - Band Width(-20dB):9MHz
 - Number of Elements: 256
 - Convex Radius: 11.7 mm
 - Volume Sweep Radius: 11.7 mm
 - FOV: 151° (B), 151*120° (Volume scan)
 - Foot Print: 29,95 (B) x 27,68 (V) mm
 - Doppler Transmission
 - Low Freq.: 6,67 MHz
 - Mid Freq.: 7,15 MHz
 - High Freq.: 8,34 MHz
 - Harmonic Transm. Frequency:
 - Resolution: 5,89 MHz
 - Normal: 5,89 MHz
 - Penetration: 5,56 MHz
 - Biopsy Guide Available: PEC63, Single-Angle, Reusable
- RSP6-16-D Wide Band Linear Volume Probe
 - Applications: Small Parts, Peripherals. Vascular, Pediatrics, Ortho
 - Maximum Band Width (-20dB): 5.6 – 18.4 MHz
 - Number of Elements: 192
 - Volume Sweep Radius: 33 mm
 - FOV: 37.4 mm (B); 37.4 mm * 29° (Volume scan)
 - Foot Print: 38.4 x 44.5 mm
 - Doppler Transm. Frequency.:
 - Low Freq.: 6,26 MHz
 - Mid Freq.: 7,15 MHz
 - High Freq.: 8,34 MHz
 - Harmonic Frequency:
 - Resolution: 6,67 MHz
 - Normal: 6,26 MHz
 - Penetration: 5,89 MHz
 - Biopsy Guide Available: PEC75, Single-Angle, Reusable

- USB (6x external, 5x internal for Dongles, etc)
- Parallel Port
- Remote BW Printer via USB
- Remote Color Printer via USB
- Remote VCR (RS232) via USB
- External microphone via USB

Safety Conformance

The Voluson E8 is:

- Listed to UL 60601-1 by a Nationally Recognized Test Lab
- Certified to CSA 22.2, 60601.1 by an SCC accredited Test Lab
- CB-Test report by National Certification Body
- CE Marked to Council Directive 93/42/EEC on Medical Devices
- Conforms to the following standards for safety:
 - EN 60601-1 Electrical medical equipment
 - EN 60601-1-1 Electrical medical equipment
 - EN 60601-1-2 Electromagnetic compatibility
 - EN 60601-1-4 Programmable medical systems
 - EN 60601-2-37 Particular requirements for the safety of ultrasound medical diagnostic and monitoring equipment
 - IEC 601157 Declaration of acoustic output
 - ISO 10993 Biological evaluation of medical devices
 - NEMA UD3 Acoustic output display (MI, TIS, TIB, TIC)
 - WEEE (Waste Electrical and Electronic Equipment)

External Inputs and Outputs

Connectivity on rear panel (direct access)

- VGA Out
- Footswitch via USB
- Network (RJ45)
- USB (2x)
- RS 232 (Optional, USB to RS232 converter)

Connectivity behind rear panel (access after opening):

- Video Out
 - RGB
 - S-Video (VTR)
- Video In:
 - S-Video
- Audio Out
 - Left/right
- Audio In
 - Left/right